



Resolution settings to install Codexpert™ on Windows 8.1 Pro (Parallels 9 / Mac OSX Mavericks)

[provided by Trish Levy, 20/02/14]


1. Shut down Window 8.1 Pro



2. Open Parallels [ensuring not to open Windows 8.1 Pro] > Virtual Machine > Configure




3. Go to the General tab and select Windows 8 from the dropdown box under the “Name:” field.



[NOTE]: This will change the features able to be selected in “Video”.

4. Go to the Hardware tab and select “Video”. Then, select “Scaled”.



[NOTE]: This will affect picture quality and will likely revert back to “Best for Retina” once installed.